# Conduct an Inventory

- I AND

## OVERVIEW

- Identify minor injury kit components (IFAK)
- Identify trauma kit components (IFAK)
- Identify primary inventory (CLS)
- Identify primary inventory (CAS)
- Identify safety precautions associated with components

## LEARNING OBJECTIVES

Please Read Your

**Terminal Learning Objectives** 

And

**Enabling Learning Objectives** 



## Individual First Aid Kit (IFAK)

- Worn by all combatants
- Intended for minor injuries
- Utilize casualty's IFAK to treat THEIR injuries
  - Inventory:
    - Adhesive Bandages
    - Burn Dressing
    - Combat Tape
    - Water Purification Tablets
    - First Aid Ointment
    - Triangular Bandages
    - TCCC Casualty Cards

- Typically supplemented with:
  - Combat Action Tourniquet
  - Combat Gauze
  - Pressure Dressings
  - Chest wound kit



## Combat Lifesaver Bag (CLS Bag)

- Worn by all CLS trained Marines
- Utilized to supplement HM medical supplies
- Inventory:
  - Adhesive bandage
  - Sterile Burn Dressing
  - Combat Tape
  - Water purification tablets
  - Bacitracin
  - Combat Gauze
  - Compressed Gauze Roll
  - "H" Bandage (Pressure Dressing)



## Combat Lifesaver Bag (CLS Bag) cont.

- Inventory:
  - Cravat
  - Chest Wound Kit
  - CAT Tourniquet
  - NPAs
  - Needle for Needle Decompression
  - Gloves
  - SAM Splints
  - TCCC CARD







Medical Assault Pack PN: MAP- 0904 NSN: 8465-01-612-4152



- Worn by HMs
- Inventory varies by mission/scope
- Primary Inventory (this list is not all inclusive):
  - CAT Tourniquets
  - Combat Gauze
  - Dressings
  - Burn Dressings
  - Adhesive Tape

- Inventory (cont.):
  - Bandages
  - Burn Dressings
  - Mole skin
  - Sterile burn blanket
  - Chest Wound Kit
  - Needle Decompression
  - NPAs
  - BP cuff

- Inventory (cont.):
- Stethoscope
- Bacitracin
- SAM splints
- Gloves
- Cric set
- IV kit
- Rigid eye shield

- Inventory (cont.):
- TCCC CARD
- Extraglottic Airway
- Intraosseous Infusion System (Talon)
- Gloves
- APLS thermal guard system
- Sharpie

#### Safety

Corpsmen are responsible for the education of all personnel on proper use and safety requirements of all IFAK and CLS bag equipment.

This includes:

- Indications/Contraindications
- Proper Application
- Allergy identification (particularly with Betadine)
- Common errors in application



## **Conduct** an Inventory

E